

DESCRIPTION OF SERVICES:

Written Respiratory Protection Program \$250.00

1.0	Purpose
2.0	Scope and Application
3.0	Responsibilities Program Administrator Supervisors Employees Program Elements
	Program Administrator
	Supervisors
	Employees
4.0	Program Elements
	Selection Procedures
	Medical Evaluation
	Fit Testing
	Fit Testing Respirator Use
	Air Quality
	Cleaning, Maintenance, Change Schedules and Storage
	Training
5.0	Program Evaluation

6.0 Documentation and Recordkeeping

On-line Medical Evaluations: \$25.00 each

- Completed on-line accessible from any internet connection
- HIPPA compliant
- One log on specific to your company
- · Results are real-time reported to the program administrator via email
- 100% access to your data on-line 24/7
- Includes phone consultation with our physician if initial review is deferred
- No set up fee
- Invoiced monthly based on usage

On-line Train-the-Trainer for Qualitative Fit Testing \$350.00/session up to 3 participants:

- Each session includes
 - o Power Point
 - o Lecture
 - o Detailed hands-on approach
 - Fit Testing/Training
 - Certificate of Completion

On-line Training:

- Instructions on Uses and Limitations
- Inspection
- Donning and Doffing
- User Seal Check
- Maintenance and Storage

Qualitative Respirator Fit Test Kit \$305.00

Technician On-site To Perform Qualitative Fit Testing (QLFT): Pricing based on volume and location

• Up to 96 tests/day. Sessions begin every 30 minutes

A QLFT is a pass/fail fit test that assesses the adequacy of respirator fit based on the individual's response to a test agent. A QLFT may be used for all tight fitting negative pressure facepieces that will be used in atmospheres up to 10 times the OSHA PEL. As part of conducting the fit test, technicians will conduct a user seal check on each respirator that is fit tested. QLFT is conducted using our fit fest kits and the test agent Saccharin Solution Aerosol or Bitrex[™] Solution Aerosol.

All employees to be fit tested must be clean-shaven. Each employee must be in possession of, and fit tested on, his or her own respirator at time of testing. Surrogates, or back-up respirators, must be provided by the client and on hand and available for use at the customer's premises on the day of testing.

Training:

- Instructions on Uses and Limitations
- Inspection
- Donning and Doffing
- User Seal Check
- Maintenance and Storage

Reports Provided:

Service	Report	Report Format	Timeline
Qualitative Fit Testing	Employer Summary Report indicating which employees passed or failed fit testing Employee Certification Card	Paper Copies left on-site by Technician	Day of Testing